

Configuration BD FACSVerser Rangueil

Blue Laser 488nm

AF 488 FITC GFP
Bodipy Calcein CFSE
Live Dead Green
TO-PRO1 TO-TO1 Cy2
CyQUANT DNA eGFP
eYFP Fluo3 Fluo4
PKH 67 Sytox Green
Mitotracker Green Dioc-6

PerCP
PerCP-Cy5.5
PE-Cy5
PE-AF647
7-AAD
(! Spectre large)

665 LP
700/54

507 LP
527/32

Blank

752 LP
783/56

PE-Cy7

560 LP
586/42

PE

PKH26 PE-Texas Red
SNARF 1(pH6) PE-CF 594
PE-AF 610
IP (! Spectre large)

488/15

SSC

